

Tutorial (2)

Held before FIT beam test @ PS (06.2016)


30.05.2016, Vidyo

Maciej Slupecki
University of Jyväskylä, Finland




UNIVERSITY OF JYVÄSKYLÄ

Part 2 – outline



- DAQ scheme
- Data collection
 - DAQ program, binary file format
 - Examples of raw waveforms
 - Which waveform qualities are interesting?
- Data analysis chain
 - Raw waveform analysis – existing solution
 - DigDataProcessing
 - Usage instructions
 - Format of output root file
 - Root file analysis
 - TBrowser (exercise)
 - Simple script: genericTreeReader (exercise)
 - More advanced program: fitTest-2016-06 (work-in-progress)
- DigDataProcessing modification

Part 2 – outline





- DAQ scheme
- Data collection
 - DAQ program, binary file format
 - Examples of raw waveforms
 - Which waveform qualities are interesting?
- Data analysis chain
 - Raw waveform analysis – existing solution
 - DigDataProcessing
 - Usage instructions
 - Format of output root file
 - Root file analysis
 - TBrowse (exercise)
 - Simple script: genericTreeReader (exercise)
 - More advanced program: fitTest-2016-06 (work-in-progress)
- DigDataProcessing modification



Part 2 – outline

- DAQ scheme
- Data collection
 - DAQ program, binary file format
 - Examples of raw waveforms
 - Which waveform qualities are interesting?
- Data analysis chain
 - Raw waveform analysis – existing solution
 - DigDataProcessing
 - Usage instructions
 - Format of output root file
 - Root file analysis
 - TBrowser (exercise)
 - Simple script: genericTreeReader (exercise)
 - More advanced program: fitTest-2016-06 (work-in-progress)
- DigDataProcessing modification

DAQ scheme



DAQ scheme




Name	T0-C	T0-D	V0+	XY	6x6 det	2MCP*	1MCP*
Detector type	PMT	PMT	PMT / MCP / SiPM	Custom	Custom	MCP	MCP
Number of signals	1	1	1-16?	2	1	4-5	4-5
Triggering?	yes	yes	no	no	(yes, if present)	no	no

*any of the following is possible: 1MCP, 2MCP, 4MCP, 5MCP


Digitizer & DAQ program



- Digitizer: CAEN DT5742
 - <http://www.caen.it/jsp/Template2/CaenProd.jsp?parent=14&idmod=651>
- DAQ program: WaveDump
 - <http://www.caen.it/jsp/Template2/CaenProd.jsp?idmod=692&parent=38>
 - Needs a configuration file (which should be prepared by experts)



WaveDump usage



[space]	Display help (full command list)
q	quit WaveDump
s	start data acquisition*
W [shift+w]	enable continuous writing to disk
p	plot a single event
P [shift+p]	enable continuous event plotting
0..7	enable or disable plotting waveform from a given channel
g	switch to the next group of 8 channels**

*** The data is not being saved to the disk at this point**


** The plotter is able to plot 8 channels at maximum, so to plot the next group of 8 channels the 'g' switch is introduced

The algorithm to start DAQ is following:

- Start WaveDump and press 's'
- Check from the terminal that some data is being collected and/or start plotting by pressing 'P'
- If all is ok, start writing to disk: 'W'


Data analysis chain

Recap




- Detectors
- FE electronics
- CAEN Digitizer
- DAQ program
- **Binary files**
- Waveform analysis
- Root file
- Final analysis

Binary data format




- Data from every channel goes to a separate file → acquisition directory contains the following files:
 - wave_0.dat, wave_1.dat, ..., wave_15.dat
- Each file contains just the data, there are no headers, trailers or consistency error bits, for example **wave_0.dat** has the following structure:
 - event 0:
 - starts at 0th byte of the file
 - ends at 4095th byte
 - contains **1024** float numbers (**f**)
 - each **f** is 4 bytes long and represents the voltage (measured in arbitrary units) in channel zero (coming from wave_0.dat file) at the time (200 ps * f_position), where f_position ranges from **0** to **1024**
 - The whole event (set of **1024 f**-numbers) forms a waveform (or oscillogram) with the total acquisition length of 204.8 ns and sampling rate of 5Gs/s
 - event 1:
 - starts at 4096th byte of the file
 - ends at 8191th byte
 - ...

Binary data format




- Data from every channel goes to a separate file → acquisition directory contains the following files:
 - wave_0.dat, wave_1.dat, ..., wave_15.dat
- Each file contains just the data, there are no headers, trailers or consistency error bits, for example **wave_0.dat** has the following structure:
 - event 0:
 - starts at 0th byte of the file
 - ends at 4095th byte
 - contains **1024** float numbers (**f**)
 - each **f** is 4 bytes long and represents the voltage (measured in arbitrary units) in channel zero (coming from wave_0.dat file) at the time (200 ps * f_position), where f_position ranges from **0** to **1024**
 - The whole event (set of **1024 f**-numbers) forms a waveform (or oscillogram) with the total acquisition length of 204.8 ns and sampling rate of 5Gs/s
 - event 1:
 - starts at 4096th byte of the file
 - ends at 8191th byte
 - ...

Binary data format



- Data from every channel goes to a separate file → acquisition directory contains the following files:
 - wave_0.dat, wave_1.dat, ..., wave_15.dat
- Each file contains just the data, there are no headers, trailers or consistency error bits, for example **wave_0.dat** has the following structure:
 - event 0:
 - starts at 0th byte of the file
 - ends at 4095th byte
 - contains **1024** float numbers (**f**)
 - each **f** is 4 bytes long and represents the voltage (measured in arbitrary units) in channel zero (coming from wave_0.dat file) at the time (200 ps * f_position), where f_position ranges from **0** to **1024**
 - The whole event (set of **1024 f**-numbers) forms a waveform (or oscillogram) with the total acquisition length of 204.8 ns and sampling rate of 5Gs/s
 - event 1:
 - starts at 4096th byte of the file
 - ends at 8191th byte
 - ...

Binary data format



- Data from every channel goes to a separate file → acquisition directory contains the following files:
 - wave_0.dat, wave_1.dat, ..., wave_15.dat
- Each file contains just the data, there are no headers, trailers or consistency error bits, for example **wave_0.dat** has the following structure:
 - event 0:
 - starts at 0th byte of the file
 - ends at 4095th byte
 - contains **1024** float numbers (**f**)
 - each **f** is 4 bytes long and represents the voltage (measured in arbitrary units) in channel zero (coming from wave_0.dat file) at the time (200 ps * f_position), where f_position ranges from **0** to **1024**
 - The whole event (set of **1024 f**-numbers) forms a waveform (or oscillogram) with the total acquisition length of 204.8 ns and sampling rate of 5Gs/s
 - event 1:
 - starts at 4096th byte of the file
 - ends at 8191th byte
 - ...

Wake up question



Keeping in mind that the sampling rate of the digitizer is 5 Gs/s → **200 ps.**

How can we study the detector **time** features (resolution) with **~10 ps accuracy**?



ALICE
A JOURNEY OF DISCOVERY


Raw waveforms



Raw waveforms

0: Find zero-level

1: Find amplitude






ALICE
A JOURNEY OF DISCOVERY

Raw waveforms

2: Set 10% and 90% amplitude thresholds






ALICE
A JOURNEY OF DISCOVERY

Raw waveforms

3: Fit the leading edge (LE)






ALICE
A JOURNEY OF DISCOVERY

Raw waveforms


4: Get time of LE at half-amplitude level



Raw waveforms

Issues to consider

- What if:
 - Amplitude too **big**
 - Amplitude too **small**
 - LE very **steep**
 - Shape **anomalous**
 - a smaller pre-pulse at the LE of main pulse
 - **Several** pulses



Data analysis chain

Recap

- Detectors
- FE electronics
- CAEN Digitizer
- DAQ program
- Binary files
- **Waveform analysis**
- Root file
- Final analysis

Waveform analysis

DigDataProcessing - compilation

- Copy digdataprocessing.zip from indico to your disk and unpack:
 - <https://indico.cern.ch/event/532602/attachments/1280341/1902644/digdataprocessing.zip>
- In terminal go to the directory where you just unpacked digdataprocessing software


```
cd /somewhere/digdataprocessing
```

```
make
```

```
./digdataprocessing
```

- if everything is working fine then the following should be printed

```
<I> === digdataprocessing ===  
<I>   - FIT / T10@CERN beam test software.  
<I>   - Takes as input raw binary files generated by CAEN Digitizer DT57**  
<I>   - Plots misc stuff and/or converts data to root tree  
  
<I> Command line syntax:  
<I>   -g          ---> GUI mode  
<I>   -c          ---> text mode, convert binary data to root  
<I>   -i=in_dir   ---> input directory path where run folders are placed  
<I>   -o=out_path ---> output file path where root file is saved [default: in_dir/../rootFiles/in_dir.root]  
<I>   -n=n_events ---> number of events to be processed [default => all events]  
  
<I> Example:  
<I>   ./digdataprocessing -g  
<I>   ./digdataprocessing -c -i=/path_data/
```



Waveform analysis

DigDataProcessing - compilation

- Copy digdataprocessing.zip from indico to your disk and unpack:
 - <https://indico.cern.ch/event/532602/attachments/1280341/1902644/digdataprocessing.zip>
- In terminal go to the directory where you just unpacked digdataprocessing software

```
cd /somewhere/digdataprocessing
```

```
make
```

```
./digdataprocessing
```

- if everything is working fine then the following should be printed

```
<I> === digdataprocessing ===  
<I>   - FIT / T10@CERN beam test software.  
<I>   - Takes as input raw binary files generated by CAEN Digitizer DT57**  
<I>   - Plots misc stuff and/or converts data to root tree  
  
<I> Command line syntax:  
<I>   -g          ---> GUI mode  
<I>   -c          ---> text mode, convert binary data to root  
<I>   -i=in_dir   ---> input directory path where run folders are placed  
<I>   -o=out_path ---> output file path where root file is saved [default: in_dir/../rootFiles/in_dir.root]  
<I>   -n=n_events ---> number of events to be processed [default => all events]  
  
<I> Example:  
<I>   ./digdataprocessing -g  
<I>   ./digdataprocessing -c -i=/path_data/
```



ALICE
A JOURNEY OF DISCOVERY

Waveform analysis

DigDataProcessing – compilation / troubleshooting

- If 'make' command fails:
 - check that the root related environmental variables are set correctly:
 - 'echo \$PATH' should contain: *root/bin/*
 - 'echo \$LD_LIBRARY_PATH' should contain: *root/lib/*
 - If any of them are missing then run:
 - 'source /path_to_your_root/bin/thisroot.sh'
 - Hint: This line can be added to your ~/.bashrc file (this file is executed every time a new terminal is opened)
 - re-run make (after navigating back to the digdataprocessing directory)
- If the digdataprocessing executable is not present after running 'make' command → it means 'make' failed
- In some cases when the source code is modified, before running 'make' try 'make clean'



ALICE
A JOURNEY OF DISCOVERY


Waveform analysis

DigDataProcessing – compilation / troubleshooting

- If 'make' command fails:
 - check that the root related environmental variables are set correctly:
 - 'echo \$PATH' should contain: *root/bin/*
 - 'echo \$LD_LIBRARY_PATH' should contain: *root/lib/*
 - If any of them are missing then run:
 - 'source /path_to_your_root/bin/thisroot.sh'
 - Hint: This line can be added to your ~/.bashrc file (this file is executed every time a new terminal is opened)
 - re-run make (after navigating back to the digdataprocessing directory)
- If the digdataprocessing executable is not present after running 'make' command → it means **'make'** failed
- In some cases when the source code is modified, before running **'make'** try **'make clean'**

Waveform analysis

DigDataProcessing – usage



- Download a sample of raw data collected previously and unpack it (1.1 Gb):
 - <https://filesender.funet.fi/?vid=297c1cf5-0de9-74c8-a3c5-00004a5a92c1>
- Assuming the data has been unpacked into: /alicedata/ run:

```
./digdataprocessing -c -i=/alicedata/08102015_1758_RUN17_  
NBI1655frag22_MOD1618frag4_6gevc_collimator_min_3x3
```

- The program stops immediately and asks that user creates a directory called **rootFiles** in **/alicedata/**, which is where the output files will be stored by default. Do it and **rerun**.
- The default output location can be overridden by adding to the command line:


```
o=/path_where_the_root_files_go/filename.root
```

Note that the **full path** and **full name** of the file has to be provided with this option

- The digdataprocessing is running → converting every waveform into event-channel amplitude-time set of numbers
 - (+ some more characteristics, like total charge, time reconstructed using different methods, etc)
- After it is finished, the result is pre-analysed **root file** (2 Gb → 10 Mb)

Waveform analysis

DigDataProcessing – usage



- Download a sample of raw data collected previously and unpack it (1.1 Gb):
 - <https://filesender.funet.fi/?vid=297c1cf5-0de9-74c8-a3c5-00004a5a92c1>
- Assuming the data has been unpacked into: /alicedata/ run:

```
./digdataprocessing -c -i=/alicedata/08102015_1758_RUN17_  
NBI1655frag22_MOD1618frag4_6gevc_collimator_min_3x3
```

- The program stops immediately and asks that user creates a directory called **rootFiles** in **/alicedata/**, which is where the output files will be stored by default. Do it and **rerun**.
- The default output location can be overridden by adding to the command line:


```
o=/path_where_the_root_files_go/filename.root
```

Note that the **full path** and **full name** of the file has to be provided with this option

- The digdataprocessing is running → converting every waveform into event-channel amplitude-time set of numbers
 - (+ some more characteristics, like total charge, time reconstructed using different methods, etc)
- After it is finished, the result is pre-analysed **root file** (2 Gb → 10 Mb)

Data analysis chain


Recap



- Detectors
- FE electronics
- CAEN Digitizer
- DAQ program
- Binary files
- Waveform analysis
- **Root file**
- Final analysis

Waveform analysis

DigDataProcessing – output root file format, tree structure



- Root file (**TFile**, “08102015_1758_RUN17...root”)
 - Root tree (**TTree**, “RawDataTree”)
 - Branches: channels related to a specific detector or pixel Examples:
 - **4MCP2** \Leftrightarrow **4MCP** is the one purchased by NBI, **2** is the pixel number
 - **T0_C1, T0_C2** \Leftrightarrow Two signals coming from the same triggering PMT (but not identical*)
 - **T0_D1, T0_D2** \Leftrightarrow The other triggering PMT (see slide 5)
 - **X, Y** \Leftrightarrow 2 position signals coming from XY detector
 - Two types of branches for each channel:
 - **main_4MCP2** – usable, stable and tested
 - **4MCP2** – usable only in case of **X** and **Y** branches, otherwise treat as experimental (not fully implemented and tested, in general avoid using it)


*Digitizer has 16 input channels in 2 groups: 8 channels each. Each group has a separate internal clock (no synch between them), so to obtain accurate timing a reference (trigger) signal is provided for each group separately. So when relative timing of given channel is calculated an appropriate pair of trigger signals should be used (either C1&D1 or C2&D2) as reference.

There are two ways to check which signal was in which channel originally:

- Config file: DataConfig.txt – located in the same folder as binary data
- Logbook (link from beam test performed in 10.2015, there will be a new one created this year):
https://docs.google.com/document/d/1ntYFb_oJMFxbDO4RM72wBJIVpJs8qH_kdIdd5wCkAzM/edit?usp=sharing

Waveform analysis

DigDataProcessing – output root file format, tree structure



- Root file (**TFile**, “08102015_1758_RUN17...root”)
 - Root tree (**TTree**, “RawDataTree”)
 - Branches: channels related to a specific detector or pixel Examples:
 - **4MCP2** \Leftrightarrow **4MCP** is the one purchased by NBI, **2** is the pixel number
 - **T0_C1, T0_C2** \Leftrightarrow Two signals coming from the same triggering PMT (but not identical*)
 - **T0_D1, T0_D2** \Leftrightarrow The other triggering PMT (see slide 5)
 - **X, Y** \Leftrightarrow 2 position signals coming from XY detector
 - Two types of branches for each channel:
 - **main_4MCP2 – usable, stable and tested**
 - **4MCP2** – usable only in case of **X** and **Y** branches, otherwise treat as experimental (not fully implemented and tested, in general avoid using it)


*Digitizer has 16 input channels in 2 groups: 8 channels each. Each group has a separate internal clock (no synch between them), so to obtain accurate timing a reference (trigger) signal is provided for each group separately. So when relative timing of given channel is calculated an appropriate pair of trigger signals should be used (either C1&D1 or C2&D2) as reference.

There are two ways to check which signal was in which channel originally:

- Config file: DataConfig.txt – located in the same folder as binary data
- Logbook (link from beam test performed in 10.2015, there will be a new one created this year):
https://docs.google.com/document/d/1ntYFb_oJMFxbDO4RM72wBJIVpJs8qH_kdIdd5wCkAzM/edit?usp=sharing

Waveform analysis

DigDataProcessing – output root file format, tree structure



- Root file (**TFile**, “08102015_1758_RUN17...root”)
 - Root tree (**TTree**, “RawDataTree”)
 - Branches: channels related to a specific detector or pixel Examples:
 - **4MCP2** \Leftrightarrow **4MCP** is the one purchased by NBI, **2** is the pixel number
 - **T0_C1, T0_C2** \Leftrightarrow Two signals coming from the same triggering PMT (but not identical*)
 - **T0_D1, T0_D2** \Leftrightarrow The other triggering PMT (see slide 5)
 - **X, Y** \Leftrightarrow 2 position signals coming from XY detector
 - Two types of branches for each channel:
 - **main_4MCP2 – usable, stable and tested**
 - **4MCP2** – usable only in case of **X** and **Y** branches, otherwise treat as experimental (not fully implemented and tested, in general avoid using it)


*Digitizer has 16 input channels in 2 groups: 8 channels each. Each group has a separate internal clock (no synch between them), so to obtain accurate timing a reference (trigger) signal is provided for each group separately. So when relative timing of given channel is calculated an appropriate pair of trigger signals should be used (either C1&D1 or C2&D2) as reference.

There are two ways to check which signal was in which channel originally:

- Config file: DataConfig.txt – located in the same folder as binary data
- Logbook (link from beam test performed in 10.2015, there will be a new one created this year):
https://docs.google.com/document/d/1ntYFb_oJMFxbDO4RM72wBJIVpJs8qH_kdIdd5wCkAzM/e/dit?usp=sharing

Waveform analysis


DigDataProcessing – output root file format, tree structure



- **Leaves** inside branches
 - each leaf corresponds to a waveform property extracted from the binary data
- inside branch **main_5MCP2**:
 - **Amplitude** – pulse amplitude (slide 15)
 - Time – default method of calculating time (?)
 - **Time_pol1** – leading edge fitted with polynomial of the 1st order (linear fitting) → most accurate
 - Time_pol3 – fitted with polynomial of the 3rd order → less reliable
 - ZeroLevel – (top of slide 15)
 - ZeroLevelRMS – quality check
 - time_peak – for internal use
 - time_begin – for internal use
 - time_end – for internal use
 - time_front_end – for internal use
 - time_back_begin – for internal use
- inside branch **X**:
 - Level – the level of the pulse (~amplitude of the flat part, see slide 14) corresponding to the stripe(s) hit
 - Xch – x-position in channels (1 stripe → 1 channel): 0-5
 - Xmm – x-position in mm

Waveform analysis


DigDataProcessing – output root file format, tree structure



- **Leaves** inside branches
 - each leaf corresponds to a waveform property extracted from the binary data
- inside branch **main_5MCP2**:
 - **Amplitude** – pulse amplitude (slide 15)
 - Time – default method of calculating time (?)
 - **Time_pol1** – leading edge fitted with polynomial of the 1st order (linear fitting) → most accurate
 - Time_pol3 – fitted with polynomial of the 3rd order → less reliable
 - ZeroLevel – (top of slide 15)
 - ZeroLevelRMS – quality check
 - time_peak – for internal use
 - time_begin – for internal use
 - time_end – for internal use
 - time_front_end – for internal use
 - time_back_begin – for internal use
- inside branch **X**:
 - Level – the level of the pulse (~amplitude of the flat part, see slide 14) corresponding to the stripe(s) hit
 - Xch – x-position in channels (1 stripe → 1 channel): 0-5
 - Xmm – x-position in mm

Waveform analysis

DigDataProcessing – output root file format, tree structure



- **Leaves** inside branches
 - each leaf corresponds to a waveform property extracted from the binary data
- inside branch **main_5MCP2**:
 - **Amplitude** – pulse amplitude (slide 15)
 - Time – default method of calculating time (?)
 - **Time_pol1** – leading edge fitted with polynomial of the 1st order (linear fitting) → most accurate
 - Time_pol3 – fitted with polynomial of the 3rd order → less reliable
 - ZeroLevel – (top of slide 15)
 - ZeroLevelRMS – quality check
 - time_peak – for internal use
 - time_begin – for internal use
 - time_end – for internal use
 - time_front_end – for internal use
 - time_back_begin – for internal use
- inside branch **X**:
 - Level – the level of the pulse (~amplitude of the flat part, see slide 14) corresponding to the stripe(s) hit
 - Xch – x-position in channels (1 stripe → 1 channel): 0-5
 - Xmm – x-position in mm

Root file analysis

TBrowser



- Quick browsing of a root file.

In terminal type:

```
root
```

```
new TBrowser
```

- Navigate through the disk using the left panel and double click the root file to open it (right)
- Then open the RawDataTree to see branch structure (down)



Root file analysis

TBrowser



- Quick browsing of a root file.


In terminal type:

```
root
```

```
new TBrowser
```

- Navigate through the disk using the left panel and double click the root file to open it (right)
- Then open the RawDataTree to see branch structure (down)





Root file analysis

TBrowser



ALICE

A JOURNEY OF DISCOVERY

Browser

File

Edit

View

Options

Tools

Help

Files

Draw Option:

alicedata

08102015_1758_RUN17_NBI1655frag22_MOD1

08102015_2230_RUN19_NBI1655frag22_MOD1

rootFiles

08102015_1758_RUN17_NBI1655frag22_M

RawDataTree;1

main_5MCP1

5MCP1

main_5MCP2

Amplitude

Time

Time_pol1

Time_pol3

ZeroLevel

ZeroLevelRMS

time_peak

time_begin

time_end

Filter: All Files (*.*)

Canvas_1


Editor 1

main_5MCP2.Amplitude

htemp	
Entries	19765
Mean	108.4
RMS	216.2

Command

Command (local):



Root file analysis

TBrowser



ALICE

A JOURNEY OF DISCOVERY

Browser

File

Edit

View

Options

Tools

Help

Files

Draw Option:

alicedata

08102015_1758_RUN17_NBI1655frag22_MOD1

08102015_2230_RUN19_NBI1655frag22_MOD1

rootFiles

08102015_1758_RUN17_NBI1655frag22_M

RawDataTree;1

main_5MCP1

5MCP1

main_5MCP2

Amplitude

Time

Time_pol1

Time_pol3

ZeroLevel

ZeroLevelRMS

time_peak

time_begin

time_end

Filter: All Files (*.*)

Canvas_1

Editor 1

main_5MCP2.Amplitude

htemp

Entries	19765
Mean	108.4
RMS	216.2

Command

Command (local):

Data analysis chain

Recap

- Detectors
- FE electronics
- CAEN Digitizer
- DAQ program
- Binary files
- Waveform analysis
- Root file
- **Final analysis**

Root file analysis

genericTreeReader


- TBrowser is great to take a quick peek at the data, but usually not sufficient to perform event-by-event analysis of multi-channel data
 - Solution?

Root file analysis

genericTreeReader

- TBrowser is great to take a quick peek at the data, but usually not sufficient to perform event-by-event analysis of multi-channel data
 - Solution?

Write own analysis script / program



Root file analysis

genericTreeReader – compilation and running



- A sample code: genericTreeReader. Download from indico and unpack:
 - <https://indico.cern.ch/event/532602/attachments/1280341/1901690/genericTreeReader.zip>
 - A simple heavily commented **example**, how to read a tree.
 - Experiment with **main.cpp** file and have fun! (3 exercises inside).
 - Start from **main()** function (line **153**). The whole relevant algorithm is in the **ReadAndDraw()** function.
 - In terminal go to the directory where you just unpacked genericTreeReader script

```
cd /somewhere/genericTreeReader
```

```
make
```

```
./genericTreeReader
```

- There will be similar message as in the case of **digdataprocessing** program
 - Provide the input file path to the file, which you created at slide 26, with: `-i=/path_to_root_file_to_be_read`
 - Or use the root file already included in the genericTreeReader directory and leave the '-i=' field empty:

```
./genericTreeReader -i=
```
 - Both of these root files should be equivalent (possibly even identical)
- Alternatively run as root script. The effect will be similar. In terminal type:

```
root
```

```
.L main.cpp
```

```
ReadAndDraw("/file_path/to_your/data_file.root");
```



ALICE
A JOURNEY OF DISCOVERY

Root file analysis

genericTreeReader – compilation and running

- A sample code: genericTreeReader. Download from indico and unpack:
 - <https://indico.cern.ch/event/532602/attachments/1280341/1901690/genericTreeReader.zip>
 - A simple heavily commented **example**, how to read a tree.
 - Experiment with **main.cpp** file and have fun! (3 exercises inside).
 - Start from **main()** function (line **153**). The whole relevant algorithm is in the **ReadAndDraw()** function.
 - In terminal go to the directory where you just unpacked genericTreeReader script

```
cd /somewhere/genericTreeReader
```

```
make
```

```
./genericTreeReader
```

- There will be similar message as in the case of **digdataprocessing** program
 - Provide the input file path to the file, which you created at slide 26, with: `-i=/path_to_root_file_to_be_read`
 - Or use the root file already included in the genericTreeReader directory and leave the '-i=' field empty:

```
./genericTreeReader -i=
```
 - Both of these root files should be equivalent (possibly even identical)
- Alternatively run as root script. The effect will be similar. In terminal type:

```
root
```

```
.L main.cpp
```

```
ReadAndDraw("/file_path/to_your/data_file.root");
```



ALICE
A JOURNEY OF DISCOVERY

Root file analysis

genericTreeReader – compilation and running

- A sample code: genericTreeReader. Download from indico and unpack:
 - <https://indico.cern.ch/event/532602/attachments/1280341/1901690/genericTreeReader.zip>
 - A simple heavily commented **example**, how to read a tree.
 - Experiment with **main.cpp** file and have fun! (3 exercises inside).
 - Start from **main()** function (line **153**). The whole relevant algorithm is in the **ReadAndDraw()** function.
 - In terminal go to the directory where you just unpacked genericTreeReader script

```
cd /somewhere/genericTreeReader
```

```
make
```

```
./genericTreeReader
```

- There will be similar message as in the case of **digdataprocessing** program
 - Provide the input file path to the file, which you created at slide 26, with: `-i=/path_to_root_file_to_be_read`
 - Or use the root file already included in the genericTreeReader directory and leave the '-i=' field empty:

```
./genericTreeReader -i=
```
 - Both of these root files should be equivalent (possibly even identical)
- Alternatively run as root script. The effect will be similar. In terminal type:

```
root
```

```
.L main.cpp
```

```
ReadAndDraw("/file_path/to_your/data_file.root");
```




ALICE
A JOURNEY OF DISCOVERY


Root file analysis

genericTreeReader – result


Amplitude histo



Time histo




2D Time-Amp



Root file analysis

Next steps




- More advanced examples of tree analysis can be found in digdataprocessing source code:
 - files: DigDataAnalysis.cpp and DigDataAnalysis.h
 - The features from there are not included yet as a switchable option when running the main program
 - the main developer of digdataprocessing and a contact person is:
 - Dima Finogeev (dmitry-finogeev@yandex.ru)
- Extend the genericTreeReader example
 - or just use it to read and select data, save it to a text file and use origin or excel (not recommended) to plot it

Advanced part: digDataProcessing modification

- Adjustments to meet the specifics of V0:
 - Will start this week (Dima)
 - It is a more complex task, but if there are people willing to contribute, we can arrange sth ...

Summary



- Detectors
- FE electronics
- CAEN Digitizer
- DAQ program
- Binary files
- Waveform analysis
- Root file
- Final analysis
- no details
- no details
- technical specs
- how to run
- structure
- examples, algorithm
- structure, reading
- example programs

And finally some dates ...

... and other info

- The beam time: 15.06 – 22.06
- I will be at CERN: 07.06 – 23.06
 - You are welcome to contact me at **maciej.slupecki@gmail.com** in case of any problems or questions, including:
 - problems with compilation or understanding any of the examples provided here
 - practical issues
 - Don't forget your power adapters.

